



Python module

to process FISH ([Fluorescence In Situ Hybridization](#)) images

Luis Aguilera
Joshua Cook


Motivation

- Extracting quantitative data from FISH images.



- Multi-step process. (Segmentation, spot detection, quantification).
 - Requires a skilled user.
 - Labor intensive.
 - Multiple steps require user input (to select thresholds).
 - User-to-user variability.
- Large and complex datasets that are:
 - Difficult to transfer data between computers (local/remote).
 - Difficult to interpret.
 - Difficult to apply (new/retrospective) changes.

Code architecture



Easy to use

- Experimental parameters

```
# Path to folder with images
data_dir # path to a folder with images.

# Cell Segmentation
channels_with_cytosol = [1,2] # list or int indicating the channels where the cytosol is detectable
channels_with_nucleus = 0 # list or int indicating the channels where the nucleus is detectable
channels_with_FISH = [1,2] # list or int with the channels with FISH spots.
diamter_nucleus = 200 # approximate nucleus size in pixels
diameter_cytosol = 250 # approximate cytosol size in pixels

# Parameters for spot detection
voxel_size_z = 500 # Microscope conversion px to nanometers in the z axis.
voxel_size_yx = 103 # Microscope conversion px to nanometers in the xy axis.
psf_z = 600 # PSF emitted by a [rna] spot in the z plan, in nanometers.
psf_yx = 150 # PSF emitted by a [rna] spot in the yx plan, in nanometers.

# Cluster Detection
minimum_spots_cluster = 7 # The number of intensities for a point to be considered as a cluster.
```

- Running the pipeline

```
dataframe_FISH = PipelineFISH(data_dir, channels_with_cytosol, channels_with_nucleus,
channels_with_FISH, diamter_nucleus, diameter_cytosol, voxel_size_z, voxel_size_yx, psf_z,
psf_yx, minimum_spots_cluster).run()
```


Easy to interpret

- Data-frame with all information (Pandas)


	image_id	cell_id	spot_id	nucleus_y	nucleus_x	nuc_area_px	cyto_area_px	cell_area_px	z	y	x	is_nuc	is_cluster	cluster_size	spot_type	is_cell_fragmented
0	0	0	0	260	587	10427	13213	22623	5	262	615	1	0	0	0	1
1	0	0	1	260	587	10427	13213	22623	5	296	607	1	0	0	0	1
2	0	0	2	260	587	10427	13213	22623	6	210	582	1	0	0	0	1
3	0	0	3	260	587	10427	13213	22623	7	273	577	1	0	0	0	1
4	0	0	4	260	587	10427	13213	22623	7	286	602	1	0	0	0	1

- Metadata for reproducibility (text file)

```
#####
AUTHOR INFORMATION
| Author: luisaguilera
| Created at: 10 Nov 2021
| Time: 17:4
| Operative System: darwin
| Hostname: Luis.local
#####
PARAMETERS USED
| channels_with_cytosol: [1, 2]
| channels_with_nucleus: 0
| channels_with_FISH: [1, 2]
| diameter_nucleus: 200
| diameter_cytosol: 250
| voxel_size_z: 500
| voxel_size_yx: 103
| psf_z: 600
| psf_yx: 150
| minimum_spots_cluster: 7
#####
FILES AND DIRECTORIES USED
| Directory path: /Users/luisaguilera/Desktop/FISH_Processing/dataBases/GAPDH_Exon_Cy5_short
| Folder name: GAPDH_Exon_Cy5_short
| Images in directory :
|   ROI003_XY1620755767_Z00_T0_merged.tif
#####
REPRODUCIBILITY
| Platform:
|   Python: 3.6.13
| Dependancies:
|   tqdm==4.62.3
```



Connection to NAS



- Configuration file (.yaml)

```
config_test.yaml.txt — Edited
user:
  username: user_name
  password: user_password
  remote_address : munsky-nas.engr.colostate.edu
  domain: engr_dom
```



```
NASConnection(config_file, share).copy_files(remote_folder, local_folder, timeout=60)
```

- Options to interact with NAS.
 - Download images to local computer.
 - Write .tex and .csv files to NAS.

Try it in Google Colab!

<https://colab.research.google.com/drive/1CQx4e5MQ0ZsZSQgqtLzVVh53dAg4uaQj?usp=sharing>